

Cp No. (63)918-9519618 urs.bacsecretariat2020@gmail.com

		Date :_	September 17, 2025
Company/Supplier Name	:		
Address	:		

Please quote your lowest price on the item/s listed below, subject to the General Conditions on the last page, stating the shortest time of delivery and submit your sealed quotation duly signed by your representative not later than ___09/22/25_ to URS SPMO, Morong Rizal.

Chairperson, BAC

NOTE: 1. ALL ENTRIES MUST BE TYPEWRITTEN

- 2. DELIVERY PERIOD WITHIN 30 CALENDAR DAYS
- 3. WARRANTY SHALL BE FOR A PERIOD OF SIX (6) MONTHS FOR SUPPLIES & MATERIALS, ONE (1) YEAR FOR EQUIPMENT, FROM DATE OF ACCEPTANCE BY THE ENTITY

- 4. PRICE VALIDITY SHALL BE FOR A PERIOD OF 45 CALENDAR DAYS
- 5. G-EPS REGISTRATION CERTIFICATE SHALL BE ATTACHED UPON SUBMISSION OF THE QUOTATION
- 6. BIDDERS SHALL SUBMIT ORIGINAL BROCHURES SHOWING CERTIFICATIONS OF THE PRODUCT BEING OFFERED

ITEM NO.	ITEM & DESCRIPTION	QTY.	UNIT	BID PRICE	TOTAL BID PRICE
1	Turbidity Meter Specifications:	1	unit(s)		
	Light source: LED pH measurement: yes Turbidity measurement: yes Calibration interval/record: yes Cal. adjustment for selected programs: no Special methods with pH and Temp.: yes User-defined progr. via LSdata Software: 100 Data storage: 1000 Sets: available Rechargeable batteries w. universal plug: optional Methods/programs: ca. 200 Cells / Cuvettes: round: 16 mm (H:91-104 mm); 28 mm Ident number: yes Interface RS232 PC-Connection: yes Method-update via Internet: yes Firmware-update via Internet: yes LabStation for lab use incl. rechargeable battery + PC Software LSdata optional				
	LSdata Stand alone optional optional Certificates: CE/ETLus/cETL Waterproof housing IP 67 Features:				
	Measuring range comprises 0.02 -1100 NTU/FNU with automatic measuring range switching Adjustable calibration intervals and GLP-compliant documentation Stray light behavior according to pharmacopoeia 5.0 Higly precise AMCO Clear® Standards The measurement and automatic 3-point calibration runs menu-guided. The calibration takes place via an AMCO Clear® standard set (0.02-10-1000 NTU). Up to 1000 data sets including sample ID numbers can be stored and later put out GLP-compliant via the PC software LSdata. The quality assurance of the measured results is supported by calibration interval setting with documentation.				
	Additional: *With at least 2 years warranty, calibration and preventive maintenance				
	nothing follows				
	-				

	LOT ABC - Php 200,000.00		
	-		
nd Model	:	Delivery Period	:
У		Price Validity	:
	CERTIFICATION	-	Printed Name / Signatur
and that the	CERTIFICATION tify that I have personally conducted this canvass price(s) quoted is/are true & correct and the the representative of the company who submitted		Printed Name / Signatur Tel. No. / Cell phone No.
and that the signature of	tify that I have personally conducted this canvass price(s) quoted is/are true & correct and the	- - -	
and that the signature of the quotatio	tify that I have personally conducted this canvass price(s) quoted is/are true & correct and the the representative of the company who submitted	-	Tel. No. / Cell phone No

URS-AF-AS-PRO-F-2017-02

Rev. 00

Effectivity Date August 15, 2017

Modified version of RFQ form (SF-GOOD-60, May 2004)

RFQ-2025-09-0686 PR-2025-08-1524 END USER: EDUARDO D. GATCHALIAN JR. Note: Mandatory documents
*Delivery Receipt w/ ATP Number
*Service/Billing/Sales Invoice
*Official/Collection Receipt
*Philgeps registered atleast RED membership
(Above 50k)
*Mayor's Permit
*Tax Clearance/ BIR 2303
*DTI Certificate